



VIZOR

Our company

Over 40 years developing and manufacturing contact lens care solutions and eye health products.



Our Portfolio

We offer the highest quality and compliant products available within the OTC Eye Care industry.

Presence in 90+ countries with a consolidated global network

90+ countries

In house
R&D, QA, and
Regulatory Affairs

100% independent laboratory, capable of developing complete products end-to-end: R&D, Regulatory Affairs, manufacturing, and packaging

11.000+m²
of production facilities

300+
employees

50+M
bottles
manufactured
per year

100%
renewable
electricity across
our facilities



Manufactured under applicable quality and regulatory requirements.



EU-MDR

FDA



Contact Lens Care Solutions

PEROXIDES

Ever Clean Plus	6
Ever Clean Scleral	8

MULTIPURPOSES

Regard	12
Alvera	14
Unica Sensitive	16
All Clean Soft	18

GPs

Multi	20
Cleaner	22
Conditioner	23

SALINES

Isosol	24
Saline	26

Eye Care Solutions

CLEANING

Lacrifresh Cleaning	32
---------------------	----

EYE FATIGUE

Lacrifresh Comfort	34
--------------------	----

EYE IRRITATION

Lacrifresh Aloe	36
-----------------	----

DRY EYE

Lacrifresh Moisture 0.10%	38
Lacrifresh Ocu-dry 0.20%	40
Lacrifresh Ocu-dry 0.30%	42
Ozean	44



Contact lens care category

Peroxide systems

Clean, remove deposits, disinfect, wet and neutralize with a tablet. Without the need to rub & rinse. Preservative free.

EVER CLEAN PLUS
EVER CLEAN SCLERAL

Multipurpose solutions

Clean, remove deposits, disinfect, wet, rinse & preserve contact lenses when they are not used.

REGARD
ALVERA
UNICA SENSITIVE
ALL CLEAN SOFT

GPs

Solutions for maintaining RGP lenses.

MULTI
CLEANER
CONDITIONER

Salines

Keep contact lenses free of irritating substances or particles.

ISOSOL
SALINE

*The true strength of a company lies in the **excellence of its products** and its relentless **drive for innovation.***



350 ml + 45 t
2 x 350 ml + 90 t
225 ml + 30 t
60 ml + 8 t



3% Hydrogen Peroxide system
with bi-layer tablet.
Eliminates Lipids & Proteins.

PRESERVATIVE
free

Benefits

- Bottle: Disinfection
- Tablet:

LIPID CLEANING
DISINFECTION



PROTEIN CLEANING
WETTING



To allow a prolonged release of catalase: Ensures disinfection and neutralisation.

Key differentiators

- **Lipase & Subtilisin: enzymes that eliminate lipid & protein deposits.**
- No need to rub & rinse.
- Preservative free.
- Complete cycle in 2 hours.



REMOVAL
OF PROTEINS
& LIPIDS



NO NEED
TO RUB
& RINSE



DISINFECTION

Users

- Demanding high efficiency in disinfection and cleaning.
- Demanding preservative-free solutions.

DISINFECTANT

3% H₂O₂

LUBRICANT

Povidone

TABLET

Enzymes (eliminate Lipids & Proteins)
+ Povidone (wetting) + Catalase (neutralising)
+ Chlorophyllin & Riboflavin

CL

Soft & RGP



300 ml + 30 t
80 ml + 8 t

EVER
clean
SCLERAL

3% Hydrogen Peroxide system
with bi-layer tablet.
Tablet with enzyme.
Effective protein elimination.

PRESERVATIVE
free



Benefits

- Bottle: Disinfection
- Tablet:

CLEANING
DISINFECTION

PROTEIN ELIMINATION
WETTING

To allow a prolonged release of catalase: Ensures disinfection and neutralisation.

Key differentiators

- **Subtilisin:** enzyme that **eliminates protein deposits.**
- No need to rub & rinse.
- Preservative free.
- Complete cycle in 2 hours.



REMOVAL
OF PROTEINS



NO NEED
TO RUB
& RINSE



DISINFECTION

Users

- Of scleral lenses with large diameter and sagittal height.
- Demanding high efficiency in disinfection and cleaning.
- Demanding preservative-free solutions.

DISINFECTANT

3% H₂O₂

LUBRICANT

Povidone

TABLET

Enzyme (eliminates Proteins) + Povidone (wetting)
+ Catalase (neutralising) + Chlorophyllin & Riboflavin

CL

Scleral



PRESERVATIVE
free



PRESERVATIVE
free



1. UPPER LAYER

With enzymes that remove lipids, proteins and other deposits adhering to contact lenses.

2. GREEN LAYER

A specific agent provides a prolonged release of catalase and controls peroxide neutralisation.

Povidone: Lubricating agent for added comfort.

Chlorophyllin and Riboflavin: Natural colour indicators.



1. UPPER LAYER

With enzyme that removes proteins and other deposits adhering to contact lenses.


2. GREEN LAYER

A specific agent provides a prolonged release of catalase and controls peroxide neutralisation.


With Povidone lubricant for enhanced user comfort.

With Riboflavin and Chlorophyllin: Natural colour indicators.


start



PLACEMENT
OF THE TABLET



DISINFECTION
EFFERVESCENCE



NEUTRALISATION
PROCESS COMPLETED




END
CONTACT LENSES READY


start

The white basket (R) is for the right eye contact lens. The blue basket contains a central space where the tablet is added.


Baskets for lenses up to 23 mm diameter and 9 mm sagittal height.



PLACEMENT
OF THE TABLET



DISINFECTION
ACTIVE PROCESS



NEUTRALISATION
PROCESS COMPLETED

CLC SOFT

multipurpose



355 ml / 100 ml

Regard®

Multipurpose solution with Oxychlorite®
Preservative free on the eye
& Active disinfectant.



7 in 1
FORMULA



SENSITIVE
EYES



DRY
EYE

Benefits

- Cleaning, lipid and protein elimination, disinfection, rinsing, moisturising, lubricating and storing.

Key differentiators

- OxyChlorite® molecule: Ensures effectiveness and safety in the lenscase.
- Preservatives when inserting CL on the eye.
- Enhanced comfort and hygiene with CLs.

Users

- With sensitive & dry eyes.
- Young people who will be using CLs for many years to come.
- With allergic reactions.

PRESERVATIVE

Oxychlorite®

LUBRICANT

Hypromellose 0.15%

CHELATING

EDTA

SURFACTANT

Poloxamer

BUFFER

Boric acid


CL

Soft



350 ml / 100 ml

alvera® Multipurpose solution with aloe vera.
Protects in the most natural way.



EYE
PROTECTION



REMOVAL
OF PROTEINS
& LIPIDS

Benefits

- Cleaning, protein and lipid elimination, disinfecting, rinsing, storing, wetting and lubricating.

Key differentiators

- Improves cleaning and reduces adhesion to the lenses of all kinds of deposits thanks to the **double chelating agents**.
- **Aloe Vera & Povidone:** advanced comfort with moisturising and lubricating benefits.

Users

- Long term contact lens users.
- Who are looking for extended comfort with CLs.
- Who produce lots of deposits.

PRESERVATIVE	ACTIVE AGENT	LUBRICANT	CHELATING	SURFACTANT	BUFFER	CL	LENSCASE
Polyhexanide 0.0002%	Aloe Vera	Povidone	EDTA + Anti protein deposits	Poloxamer	Borax/Boric	Soft	Flat



350 ml / 100 ml

Unica[®] sensitive

Multipurpose solution
with Sodium hyaluronate.
For users with sensitive eyes.

Benefits


- Cleaning, protein and lipid elimination, disinfecting, rinsing, storing, wetting and lubricating.

Key differentiators

- Enhanced moisturising thanks to double wetting agent: Sodium Hyaluronate & Povidone.
- Minimum concentration of preservatives for sensitive eyes.

Users

- Those seeking long-lasting comfort with CLs.
- With sensitive eyes.



DOUBLE
WETTING



SENSITIVE
EYES

PRESERVATIVE	LUBRICANTS	CHELATING	SURFACTANT	BUFFER	CL	LENSCASE
Polyhexanide 0.0001%	Sodium Hyaluronate Povidone	EDTA	Poloxamer	Borax/Boric	Soft	Flat

multipurpose CLC SOFT



350 ml / 100 ml



Multipurpose solution with anti protein deposits agent. Prevention and removal of protein buildup.



CLEANING



DISINFECTION



WETTING

Benefits

- Cleaning, protein and lipid removal, disinfecting, rinsing, storing and wetting.

Key differentiators

- Improves cleaning and reduces deposit adhesion through a combination of a surfactant and a double chelating agent.
- Prevents the accumulation of protein deposits on the lens surface and provides better lens maintenance.
- The chelating agents act to strengthen the antimicrobial efficacy.

Users

- Those who produce more deposits on their CLs.

PRESERVATIVE	LUBRICANT	CHELATING	SURFACTANT	BUFFER	CL	LENSCASE
Polyhexanide 0.0002%	Povidone	EDTA + Anti protein deposits	Poloxamer	Borax/Boric	Soft	Flat

RGP solution

CLC
RGP



240 ml / 120 ml / 60 ml



Multipurpose solution for RGP lenses.
Dual wetting agents:
increased comfort.



6 in 1
FORMULA



CLEANING



DISINFECTION



WETTING

Benefits

- Cleaning, protein and lipid removal, disinfecting, rinsing, storing and wetting.

Key differentiators

- Enhanced comfort thanks to Povidone and Hypromellose, lubricating & viscosifying agents.
- Thanks to the combined action of these two ingredients, the solution forms a coating on the contact lens, keeping the water molecules around its surface and retarding the loss of moisture.

Users

- Who require effective cleaning and comfort.

PRESERVATIVE

Polyhexanide
0.0002%

LUBRICANTS

Povidone
Hypromellose

CHELATING

EDTA

SURFACTANT

Poloxamer

BUFFER

Borax/Boric

CL

RGP

LENSCASE

Flat



Non-abrasive cleaning solution.

Benefits

- Keeps contact lenses clean.

Users

- Those who need effective cleaning and disinfection.



CLEANING



AGENT

Lauryl ether sulphate
Propylene glycol

STABILISER

Acrylic Acid

CL

RGP (Ortho-K & scleral)



Conditioning solution with dual wetting agent.

Benefits

- Disinfects, moisturises, rinses and preserves.

Users

- Those who need effective cleaning and disinfection.



WETTING



PRESERVATIVE

Polyhexanide

LUBRICANTS

Povidone
Hypromellose

CHELATING

EDTA

SURFACTANT

Poloxamer

BUFFER

Borax/Boric

CL

RGP (Ortho-K & scleral)



30 x 5 ml



Preservative free saline solution
with electrolytes.
Insertion and rinsing of scleral lenses.

**PRESERVATIVE
free**

Active Ingredients

- **Electrolytes:** Potassium, Calcium, Magnesium & Sodium.

Benefits

- Filling scleral lenses before insertion
- Improves comfort and time of use of contact lenses.
- Removal of irritating substances or particles and relief of ocular irritation.

Users

- Scleral lens users.
- Patients with ocular irritation symptoms.

ELECTROLYTES

Potassium
Calcium
Magnesium
Sodium



PRESERVATIVE

Free

OTHERS

Electrolytes
(Ca²⁺, K⁺, Mg²⁺)

BUFFER

Borax/Boric

CL

RGP (Scleral) & Soft

saline Ophthalmic solution for eye hygiene,
eye wash and contact lens rinsing.



350 ml / 30 x 5 ml

Benefits

- Rinses and moisturises.
- Preserves all types of contact lenses (format 350ml).

Users

- With symptoms of eye irritation.
- Of all types of contact lenses.

**PRESERVATIVE
free**



DAILY LENSES



EYE IRRITATION



ALL TYPES CL

PRESERVATIVE

Polyhexanide 350ml
Free : Unidose

CHELATING

EDTA

BUFFER

Borax/Boric

CL

Soft & RGP

A close-up of a human eye, with the iris and pupil reflecting a scenic landscape of green trees and a blue sky with clouds. The eye is set against a dark blue background.

Eye care category

Cleaning

Lacrifresh Cleaning

Eye
fatigue

Lacrifresh Comfort

Eye
irritation

Lacrifresh Aloe

Dry
eye

Moisture 0.10%

Ocu-dry 0.20%

Ocu-dry 0.30%

Ozean

*We bring **meaningful innovation**
to market, not just for ourselves,
but for our partners.*

Avizor drops range EYE CARE



PRESERVATIVE
free



PRESERVATIVE
free

PRESERVATIVE
free

PRESERVATIVE
free

PRESERVATIVE
free

SIZE (ml)

15 ml

.....

ACTIVE
MOLECULES

Povidone
& Poloxamer

OSMOLARITY
(mOsm/L)

300
(Isotonic)

.....

VISCOSITY (CP)

3-5 (Similar
to human tear)

.....

USEFUL
LIFETIME

60 d

CL
FRIENDLY

.....



lacrifresh
cleaning

EYE
FATIGUE

15 ml

.....

Povidone

300
(Isotonic)

.....

3-5 (Similar
to human tear)

.....

60 d

.....



lacrifresh
comfort

EYE
IRRITATION

10 ml

.....

Aloe vera
Hypromellose

Electrolytes
(Ca²⁺, K⁺, Mg²⁺)

300
(Isotonic)

.....

3-5 (Similar
to human tear)

.....

90 d

.....



lacrifresh
aloe

OCCASIONAL
DRY EYE

15 ml / 10 ml / 20 x 0.4 ml

.....

Sodium hyaluronate
0.10%

300
(Isotonic)

.....

3-5 (Similar
to human tear)

.....

60 d / 90 d / 1 single use

.....



lacrifresh
moisture

MILD-MODERATE
DRY EYE

10 ml / 20 x 0.4 ml

.....

Sodium hyaluronate
0.20%

Glycerin

Povidone

Electrolytes
(Ca²⁺, K⁺, Mg²⁺)

220
(Hypotonic)

.....

16-19 (Viscosity)

.....

90 d / 1 single use

.....



lacrifresh
ocu-dry

MODERATE-SEVERE
DRY EYE

10 ml / 20 x 0.4 ml

.....

Sodium hyaluronate
0.30%

Glycerin

Povidone

Electrolytes
(Ca²⁺, K⁺, Mg²⁺)

190
(+ Hypotonic)

.....

48-52 (+ Viscosity)

.....

90 d / 1 single use

.....



lacrifresh
ocu-dry

OCULAR DRYNESS
& EYE FATIGUE

10 ml

.....

Cycarthyn®

Hypromellose

Electrolytes
(Ca²⁺, K⁺, Mg²⁺)

220
(Hypotonic)

.....

3-5 (Similar
to human tear)

.....

90 d

.....



OZ ≈ AN®



Moisturising solution that cleans,
lubricates and moisturises
all types of CL.

PRESERVATIVE	LUBRICANT	SURFACTANT	VISCOSITY	OSMOLARITY	USEFUL LIFETIME	CL
Polyhexanide 0.0002%	Povidone	Poloxamer	3-5 cP	Isotonic (300mOsm/L)	60 days	All



Multidose 15 ml

lacrifresh[®]
comfort

Moisturising solution that instantly
relieves and refresh the eyes.

Benefits

- Wetting and lubricating the ocular surface.
- Rewetting and lubricating all types of contact lenses.

Key differentiators

- Provides immediate relief from **eye fatigue**.
- Lubricates and moisturises the ocular surface in situations of occasional dryness.

Users

- With eye fatigue.
- Wearers of all types of contact lenses.



LUBRICATING



EYE FATIGUE



ALL TYPES CL

PRESERVATIVE	LUBRICANT	CHELATING	BUFFER	USEFUL LIFETIME	CL
Polyhexanide 0.0002%	Povidone	EDTA	Borax/Boric	60 days	All



Multidose 10 ml

lacrifresh[®]
aloe

Multidose preservative-free solution
with Aloe vera.

Benefits

- Wetting and lubricating the ocular surface.
- Rewetting and lubricating all types of contact lenses.

Key differentiators

- **Reduces the symptoms of irritation** and dry eyes after one and three months of use.
- **Provides protection to epithelial cells** against eye irritation thanks to its specific content in aloe vera.
- **Increases the stability of the tear film** thanks to the combination of **Hypromellose** and aloe vera (DUAL ACTION).

Users

- With dry eye, eye irritation and eye fatigue.
- Wearers of all types of contact lenses.

PRESERVATIVE
free



LUBRICATING



DRY EYE
EYE IRRITATION
EYE FATIGUE



ALL TYPES CL

PRESERVATIVE	LUBRICANTS	VISCOSITY	OSMOLARITY	USEFUL LIFETIME	CL
Free	Aloe vera & Hypromellose	3-5 cP	Isotonic (300mOsm/L)	90 days	All



Multidose 15 ml
Multidose 10 ml
Unidose 20 x 0.4 ml

lacrifresh® moisture

Multidose solution with Sodium Hyaluronate 0.1%.
For mild dryness.

PRESERVATIVE
free

Benefits

- Wetting and lubricating the ocular surface.
Rewetting and lubricating all types of contact lenses.

Key differentiators

- Intended to **relieve** the symptoms of **mild dryness and eye fatigue**.
- Lubricates contact lenses.
- **Preservative-free formats available** (10 ml and single dose 20 x 0.4 ml).

Users

- With mild dry eye and eye fatigue.
- Users of all types of contact lenses.



LUBRICATING



DRY EYE



ALL TYPES CL

PRESERVATIVE	LUBRICANT	VISCOSITY	OSMOLARITY	USEFUL LIFETIME	CL
Polyhexanide 0.0001% Free: 10 ml & unidose	Sodium Hyaluronate 0.10%	3-5 cP	Isotonic (300mOsm/L)	60 days : 15 ml 90 days : 10 ml 1 day : unidose	All



Multidose 10 ml
Unidose 20 x 0.4 ml

lacrifresh[®]
ocu-dry
0.20%

Multidose solution with Sodium Hyaluronate 0.2%.
For mild to moderate dryness.

PRESERVATIVE
free

Benefits

- Wetting and lubricating the ocular surface.
- Rewetting and lubricating all types of contact lenses.

Key differentiators

- Relieves the symptoms of **mild to moderate dry eye**.
- Significant decrease in corneal and conjunctival staining.

Users

- With mild to moderate dry eye.
- Users of all types of contact lenses.



LUBRICATING




DRY EYE



ALL TYPES CL

PRESERVATIVE	LUBRICANTS	OTHERS	VISCOSITY	OSMOLARITY	USEFUL LIFETIME	CL
Free	Sodium Hyaluronate 0.20%, Glycerin, Povidone	Electrolytes (Ca ²⁺ , K ⁺ , Mg ²⁺)	16-19 cP	Hypotonic (220mOsm/L)	90 days : 10ml 1 day : Unidose	All



Multidose solution with Sodium Hyaluronate 0.3%.
For moderate to severe dryness.

- Wetting and lubricating the ocular surface.
Rewetting and lubricating all types of contact lenses.

- Relieves the symptoms of **moderate to severe dry eye**.
- Improves the subjective comfort of patients and reduces corneal staining.

- With moderate to severe dry eye.
- Users of all types of contact lenses.

LUBRICATING

42



ALL TYPES CL



Multidose 10 ml

OZAN®

Multidose preservative free solution.
Relief of ocular dryness
and eye fatigue.

Benefits

- Wetting and lubricating the ocular surface.
- Rewetting and lubricating all types of contact lenses.

Key differentiators

- **PREMIUM solution** of choice **for any discomfort** of dry eye or ocular fatigue.
- Promotes **regeneration of the tear film**.


Users

- With mild to moderate dry eye.
- Users in adverse environmental conditions.
- Users of all types of contact lenses.

PRESERVATIVE
free



LUBRICATING



DRY EYE / EYE FATIGUE



ALL TYPES CL

PRESERVATIVE	LUBRICANTS	OTHERS	VISCOSITY	OSMOLARITY	USEFUL LIFETIME	CL
Free	Cycarthyn®, Hypromellose	Electrolytes (Ca ²⁺ , K ⁺ , Mg ²⁺)	3-5 cP	Hypotonic (220mOsm/L)	90 days : 10ml	All

40+ Taking Care
YEARS of Vision Worldwide

AVIZOR



Avenida de la Innovación, 2, Parque Leganés Tecnológico, 28919 Leganés Madrid

T. 902 010 678 F. 900 110 031 E. avizor@avizor.es

www.avizor.com

v01/2026